

SAF-T-VAC®

Smoke and Fluid Evacuator

Item Number: ST9000S (10 each per box)

The Ideal Solution For Removing Smoke and Fluid

Electrocauterual smoke is hazardous to your health and the health of your O.R. staff. Toxins and pathogens are carried by the smoke plume.

REMOVE IT AT THE SOURCE!

The **SAF-T-VAC** is a cost-effective and convenient method for keeping the air in the operating room clean.

| Ideal Smoke Evacuation | SAF-T-VAC |
|--|-----------|
| Removes smoke and fluids continuously at the source | ✓ |
| Sleek, unobtrusive design improves visualization | ✓ |
| Inexpensive and disposable, Saves time and money | ✓ |
| Compatible with standard and extended electrocauterual scalpel | ✓ |

To Order:

Call: 866.225.5785

On-Line: www.Surgiform.com

Fax: 803.462.1714

Item #: ST9000S (10 each per box)

- No additional new equipment to purchase.
- No capital expenditure required.
- Works with existing systems and infrastructure (operating room/suite wall suction system).



SAF-T-VAC
Electrocauterual Scalpel Sold Separately

SUCTION
CANISTER

IN-LINE
FILTER
Optional/Suggested

WALL
SUCTION

Clean, Safe and Effective Removal of Smoke and Fluid

- SAF-T-VAC can be used with blade or needle electrodes.
- Removes smoke and fluid at the source.
- Unique design allows for 360° positioning of tip and evacuation tube.
- Does not impede visibility and maneuverability.
- Ease of set-up and use for surgeon and O.R. staff.



Electrocauterual Scalpel Sold Separately

Multi-Port
Suction

Electrode
Housing

Electrode
Adapter

Electrode