

Silicone Implants



Carvable Silicone Block


Surgiform offers a comprehensive selection of silicone facial implants. All implants are expertly designed and fabricated to deliver reliable cosmetic results. With one of the largest selections to choose from, you can be sure that we have the right style and size implant for your specific needs. Our implant grade silicone elastomer produces a softer, more supple implant, giving you superior performance and long-lasting results you and your patients can count on.

Ideal Silicone Implant Features

Surgiform


Highest quality implantable silicone	✓
Provided sterile	✓
Multiple chin, malar and nasal designs	✓
Over 20 year clinical history with more than 30,000 patients	✓

STANDARD CHIN




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)
SC-01	3.2	.50	1.3
SC-02	3.6	.60	1.5
SC-03	3.9	.70	1.6
SC-05	4.5	.90	1.7

CONCAVE CHIN




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)
CC-01C	3.2	.40	1.2
CC-02C	3.4	.50	1.4
CC-03C	3.9	.60	1.6
CC-04C	4.2	.70	1.7
CC-05C	4.7	.85	1.8

DEEP CONCAVE CHIN




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)
DCC-1XC	3.5	.45	1.4
DCC-2XC	3.6	.52	1.4
DCC-4XC	4.1	.70	1.6

EXTENDED CHIN




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)
EC-01	7.6	.50	1.6	3.9
EC-02	7.6	.70	1.6	4.4
EC-03	7.6	.80	1.7	4.8
EC-04	7.6	.90	1.7	4.9

EXTENDED ANATOMICAL




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)
EAC-01	6.1	.60	1.2	.27
EAC-02	6.3	.70	1.3	.27
EAC-03	6.3	.80	1.4	.27
EAC-04	6.3	.90	1.4	.27

NASAL DORSAL SERIES A




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)
ND-01A	5.80	.75	3.50	.20
ND-02A	5.80	.90	3.50	.30
ND-03A	5.80	.90	3.50	.40
ND-04A	5.80	1.10	3.50	.50
ND-05A	5.80	1.50	3.50	.60

NASAL DORSAL SERIES AC




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)
CND-1AC	5.80	1.00	3.50	.20
CND-2AC	5.80	1.00	3.50	.30
CND-3AC	5.80	1.00	3.50	.40
CND-4AC	5.80	1.00	3.50	.50
CND-5AC	5.80	1.00	3.50	.60

NASAL DORSAL SERIES B




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)
ND-07B	5.80	.80	3.50	.50
ND-08B	5.80	1.00	3.50	.60
ND-09B	5.80	1.10	3.50	.70
ND-10B	5.80	1.20	3.50	.60

NASAL DORSAL SERIES BC




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)
CND-9BC	5.80	1.00	3.50	.50
CND-10BC	5.80	1.00	3.50	.60
CND-11BC	5.80	1.00	3.50	.70

NASAL DORSAL SERIES C




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)	Dim E (cm)
ND-11C	5.80	1.00	3.50	.40	.20
ND-12C	5.80	1.00	3.50	.50	.30
ND-13C	5.80	1.00	3.50	.50	.40
ND-14C	5.80	1.00	3.50	.60	.20

NASAL DORSAL SERIES CC




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)	Dim E (cm)
CND-11CC	5.80	1.00	3.50	.40	.20
CND-16CC	5.80	1.10	3.50	.60	.50

NASAL - OFODILE SERIES I




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)
OFD1-01	3.5	1.0	0.8	0.2
OFD1-02	3.5	1.0	0.8	0.4
OFD1-03	3.5	1.0	0.8	0.6

NASAL - OFODILE SERIES II




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)	Dim E (cm)
OFD2-01	3.5	1.0	0.8	0.2	2.5
OFD2-02	3.5	1.0	0.8	0.4	2.5
OFD2-03	3.5	1.0	0.8	0.6	2.5

NASAL DORSAL SERIES COLUMELLA




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)	Dim E (cm)
DC-06D	4.70	3.00	.55	1.10	.60
DC-07D	4.80	3.00	.60	1.20	.60
DC-08D	4.80	3.10	.65	1.20	.60
DC-09D	4.80	3.00	.45	1.10	.60
DC-10D	4.70	2.90	.50	1.10	.60

NASAL CRESCENT DORSAL SERIES COLUMELLA




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)	Dim E (cm)
CDC-06D	4.70	3.00	.55	1.10	.60
CDC-08D	4.80	3.10	.65	1.20	.60
CDC-09D	4.80	3.00	.45	1.10	.60
CDC-10D	4.80	2.90	.50	1.10	.60

ASIAN NASAL IMPLANTS - FEMALE, MALE




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)	Dim E (cm)
DC-11DF	5.70	2.50	.40	.95	.30
DC-12DM	5.90	2.50	.30	1.00	.30

ASIAN NASAL IMPLANTS - STANDARD




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)	Dim D (cm)	Dim E (cm)
DC-11D	5.20	3.50	.35	.90	.27
DC-12D	5.50	3.50	.30	1.00	.30
DC-13D	5.50	3.70	.40	1.00	.37
DC-14D	5.50	3.90	.50	1.10	.35
DC-15D	5.90	3.90	.40	1.00	.40

STANDARD MALAR




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)
SM-01	2.0	2.7	.40
SM-02	2.2	3.1	.50
SM-03	2.5	3.9	.60
SM-04	2.8	4.5	.50
SM-05	3.4	5.4	.60

EXTENDED TRAPEZOIDAL MALAR




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)
ETM-01	4.6	1.8	.40
ETM-02	4.9	2.0	.50
ETM-03	5.2	2.2	.55
ETM-04	5.5	2.2	.60

ANATOMICAL MALAR



Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)
AM-01	3.0	4.2	.32
AM-02	3.4	4.6	.36
AM-03	3.6	5.0	.40
AM-04	4.0	5.4	.44

TRAPEZOIDAL MALAR




Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)
TM-01	4.1	2.1	.35
TM-02	4.5	2.2	.40
TM-03	5.0	2.4	.45
TM-04	5.4	2.8	.50

CARVABLE SILICONE BLOCKS

Surgiform also offers a complete line of carvable silicone blocks. Each block has the ability to be carved into individual custom implant shapes. We offer three sizes of blocks that are non-sterile and available for immediate shipment.

CARVABLE SILICONE BLOCKS



Item Code	Dim A (cm)	Dim B (cm)	Dim C (cm)
SB-1	8.0	5.0	1.5
SB-2	10.0	7.0	2.5
SB-3	12.0	8.0	3.5